

# ROBOOTER

R E B O O T Y O U R L I F E

## ROBOOTER X40

TECHNOLOGY MEETS  
FUNCTIONAL ELEGANCE  
A TRUE ORIGINAL



DESIGN  
AWARD  
2022



[www.roobooter.com](http://www.roobooter.com)

# ROBOOTER X40

## BUILT FOR YOU

Ergonomic & Aesthetically Pleasing Design

- Independent shock suspension system
- Memory foam seating
- Auto-folding mechanism
- Smart control with app integration
- Monitoring with fault alerts



Speed  
Settings



Cruise  
Control



Remote  
Drive




Voice  
Settings



Battery  
Monitoring



Gyroscope  
for Stability



## DESIGNED TO GO

Everyday Functionality



Tight turning radius

Omni-wheel built for  
rough terrain

Auto-folds to fit in  
small places

Removable lithium  
battery

Built-in lighting for safety

## PRODUCT PARAMETERS

Folding	Electric Folding	Battery capacity	20Ah(Optional 10ah)
Maximum speed	3.75mph	Charging Time	4-6H
Maximum obstacle clearance	1.77"	Overall dimension	(L) 40"x (W) 23"x (H) 36"
Maximum gap clearance	4"	Folding shape dimension	(L) 33"x (W) 23"x (H) 15"
Maximum climbing angle	10°	Seat wide	16.5"
Weight without battery	92lbs.	Maximum range	15miles
Battery number	1	Side Arms	Angle Adjustable
Battery weight (20ah)	7.5lbs.	Maximum Load	120kg
Specifications of front and rear wheels	Front: 10-in. / Rear 10-in.	Recommended Ambient Temperature	-10°C ~ 50°C
Maximum capacity	≤265lbs.	Waterproof Grade	IPX4