

ROBOOTER

R E B O O T Y O U R L I F E

ROBOOTER E40

A SMART WHEELCHAIR
WITH THE ULTIMATE EXPERIENCE



www.robooter.com

ROBOTEER E40

ULTIMATE ADVENTURE
DELIGHTFUL INNOVATION
YOUR ROBOTEER E40



ENDURANCE YOU DESERVE POWERFUL & AFFORDABLE

14 Mile range with 20ah battery



ELEGANCE MEETS FUNCTIONALITY DESIGNED FOR YOU



NEXT GENERATION SOFTWARE AND HARDWARE INTEGRATION

Drive modes, automatic upgrades,
and self-monitoring alerts; drive at ease.



CLASSIC WHITE



CLASSIC GREEN

PRODUCT PARAMETERS

Folding	Manual folding	Front Wheel Material	Honeycomb PU solid
Maximum speed	4.4mph	Rear Wheel Material	Pneumatic tubeless
Maximum obstacle clearance	1.6"	Maximum capacity	330lbs.
Maximum gap clearance	4"	Battery capacity	20Ah(Optional 10ah)
Maximum climbing angle	9°	Overall dimension	(L)29"x(W)24"x(H)36"
Weight without battery	47.5lbs.	Folding shape dimension	(L)27"x(w)24"x(H)17"
Battery number	1	Seat width	18"
Battery weight (20ah)	7.5lbs.	Maximum range	14miles (7miles with 10ah battery)
Specifications of front and rear wheels	Front: 8-in. / Rear 10-in.	Color choices	Green, White